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
**Please notice:**

In the text of these instructions square brackets and bold cypers = [00] refer to the parts of the bed insert.

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### PREFACE

Dear customer,  
 Thank you for the confidence you have shown in BURMEIER.  
 The vertically adjustable bed insert LIPPE 120, 140 is factory-checked for electrical safety and function and has left our factory in faultless condition.  
 Please take time to read this manual carefully.  
 This applies especially to the assembly phase but also to daily use.  
 It is also a practical reference guide. Keep the instructions in a handy place.  
 These instructions provide information for the operator and user regarding the convenient handling and safe operation.  
 We wish you and your users much success in care-giving and we are convinced that our product will make a positive contribution.  
 BURMEIER GmbH & Co. KG

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**GENERAL NOTES**

The bed insert LIPPE 120, 140 is available in various models.

**Overview LIPPE 120, 140 versions of model:**

*Model number: 30.9956.12*

*Type 241600*

*Reclining surface 200 x 120 cm*

*Model number: 30.9956.14*

*Type 241600*

*Reclining surface 200 x 140 cm*

**Before putting the bed insert into operation for the first time:**

■ Read these instructions carefully in order to avoid damage or incorrect operation.

■ Clean and disinfect the bed insert before first use.

The bed insert LIPPE 120, 140 complies with all requirements of the guidelines 93/42/EWG for medical products. It is classified as a dynamic medical product of the class I in accordance with "Medizinproduktegesetz (German abbreviation: MPG § 13, Medical Products Act)."

Observe your obligations as operator in accordance with medical devices "Medizinprodukte-Betreiberverordnung (German abbreviation: MPBetreibV, Operators of Medical Products Ordinance)" in order to ensure a dependable and safe operation of this medical device without endangering patients, users and third parties.

The bed insert LIPPE 120, 140 is checked by an independent examining institute. As with all technical, electrical devices improper handling can lead to damage and/or injury.

**This manual contains safety guidelines, which must be observed. All persons who work with the bed insert LIPPE 120, 140 must familiarize themselves with these instructions and follow the safety guidelines.**

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## DEFINITIONS OF INVOLVED PERSONS

In this manual the following persons are involved:

### Operator

Operators (e.g.: Medical Supply Houses, Dealers and representatives, health insurance companies) are natural persons or legal entities, who utilize the bed insert LIPPE 120, 140 or authorize its use. The operator is responsible for instructing the user in the proper use of the unit.

### User

Users are persons who as a result of training and experience are authorized to assemble, perform maintenance and/or operate the the bed insert LIPPE 120, 140. Furthermore, the user can recognize and avoid possible dangers, and judge the clinical condition of the patient.

### Patient

In this manual a patient is an infirm or handicapped person who reclines in a bed with this bed insert.

### Specialist staff

Specialist staff are employees of the operator who are authorized as a result of their education and training to deliver, assemble, dismantle and transport the bed insert. Furthermore, they are familiar with the regulations for cleaning and disinfecting the unit.

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## SAFETY INSTRUCTIONS

The bed insert LIPPE 120, 140 represents the latest technological development at the time of delivery and is checked by an independent examining institute.

The bed insert LIPPE 120, 140 must not be modified.

■ Do not use the bed insert LIPPE 120, 140 if it is damaged.

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## MEANING OF THE SAFETY SYMBOLS

In these instructions the following safety symbols are used:

### Warning about injuries to persons

Dangerous voltage. Life Threatening.

General danger. Injury or life hazard.

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### Warning about property damages

Possibility of damage to drive, material or environment.

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### Other symbols

Useful tip. For easier operation or better understanding of the unit.

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**The safety symbol used does not replace the text of the safety note. Therefore read the safety instructions and follow them exactly!**

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### SAFETY INSTRUCTIONS FOR OPERATORS

- Remember your responsibilities in accordance with the MPBetreibV in order to ensure the safe operation of this unit without endangering patients, users and third parties.
- Before first use, instruct every operator in the safe operation of the bed insert in accordance with these safety instructions, which must be supplied along with the unit.
- Call attention to the dangers of improper use of the unit, especially in regards to the electrical drives and sideguards.
- The bed insert should only be operated by properly trained personnel.
- Make sure that representatives are also familiar with the operation of the bed insert.
- During long term use it is recommended that the bed insert be inspected for function and visible damage at regular intervals (Recommendation: yearly) (see chapter 6).
- Connect the bed insert directly to an electrical outlet. Do not use extension cords and/or multiple outlet extension cords.
- When attaching other equipment (e.g. compressors for positioning systems) make sure that they are mounted and can function in a secure and safe manner.
- Do not put multiple outlet electrical extensions under the bed insert. Leaking liquids can cause a fire hazard.

#### ■ Pay special attention to:

Secure locating of all wiring, cables, tubes, etc.

▶ See page 10: *Intended use.*

▶ See page 18: *Safety note for the electric drive system*

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For further information please contact the manufacturer of the accessories or BURMEIER directly.

■ Make sure that your personnel follow the safety instructions.

#### SAFETY INSTRUCTIONS FOR USERS

■ Let the operator instruct you in the safe use of the bed insert.

■ Before each use make sure that the bed insert is in proper and faultless condition.

■ Take care that no obstacles, such as furniture or slanted ceilings interfere with the adjustment functions.

■ Pay attention that when using additional electrical components, such as patient lifts, reading lights or compressors for positioning systems that their electrical cords do not get entangled or damaged by the moving parts of the bed insert.

■ Make sure that when attaching other appliances (e.g. compressors for positioning systems) that a safe mounting and function is guaranteed.

■ Do not put multiple electrical outlets under the bed insert. Leaking liquids can be a fire hazard.

■ Pay special attention to securely locating all wiring, cables, tubes, etc.

For further information please contact the manufacturer of the accessories or BURMEIER directly.

■ Take the bed insert out of operation if damage or a malfunction is suspected:

- ▶ Immediately unplug the electrical plug from the wall outlet;
- ▶ Clearly mark the bed insert as " DEFECTIVE " ;
- ▶ Promptly notify the appropriate operator.

A checklist for assessing the proper condition is located in Chapter 6.

■ Position the main cable in such a way that during normal operation of the bed insert the cable will not be stretched, driven over or interfere with moving parts of the bed insert.

■ Check the main cable regularly for damage (abrasions, exposed wires, kinks, pressure marks, etc.) in particular:

- ▶ after every larger mechanical strain e.g.: crushing with the base frame or rolling over the power cord with an equipment cart.
- ▶ after relocating/moving and before plugging in the main cable.
- ▶ during prolonged operating by the user.

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- Make sure that the low voltage connection cable of the power adapter is not damaged by the base frame or an equipment cart.
- Do not use the area under the bed insert as a "parking spot" for any utensils.
- Lower the reclining surface in a lower position when leaving the patient unattended. This reduces the risk of injury to the patient getting in and out of the bed.
- Keep the hand switch safe from accidentally falling during non-use (Hang it on the hook). Take care that the cable is not damaged by the moving parts of the bed insert.
- To protect the patient and especially children from unintentionally adjusting the hand switch, place it beyond their reach (e.g. at the foot end) or lock the hand switch when:
  - ▶ the patient is not able to control the bed insert safely or is unable to get out of dangerous positions without help;
  - ▶ the patient could be endangered, because of unintentional adjustments to the drive motor;
  - ▶ unsupervised children are in the room with the bed insert.
- Adjustment functions can only be performed by a trained person or in the presence of a trained person.
- The failure of a lifting motor can result in the bed moving into the Trendelenburg position (head end lowered), which could endanger the patient. If this occurs, the bed must be returned to a horizontal position. Lock the lifting function and inform the operator immediately.

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2,5 cm minimum tolerance between reclining surface frame and bed box

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- The reclining surface frame must not dip below the bed box in which it was installed. There is a danger that limbs could be squeezed! When lowering leave 2.5 cm tolerance. Note the safety instructions in chapters:
  - ▶ Safety guidelines for electrical drive systems (see page 18)
  - ▶ Cleaning and disinfecting (see page 26)
  - ▶ Safety guidelines for exchanging electrical components (see page 40)
 as well as safety guidelines.

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
## PRODUCT DESCRIPTION

### NORMAL OPERATING PROCEDURE

■ The bed insert LIPPE 120, 140, hereafter called bed insert, serves as a comfortable solution for positioning and easing the care of infirm, frail people in senior or nursing homes. Additionally it was developed as a comfortable solution for in-home care of infirm, handicapped or frail people. It will enhance such this care.

■ The Lippe 120, 140 bed insert is designed for fitting to existing beds or bed frames and must only be used in conjunction with these.

■ The bed insert is designed for use in hospitals only in rooms of application group 0.

■ This bed insert can be used for care under the instructions of a physician and can serve in diagnosing, treatment and observation of a patient. It is therefore equipped with a locking mechanism for the electrical adjustment functions.

■ This bed insert has no special connector for potential equalization. Consider this before combining with additional electrically operated (medical) devices. Further references about additionally applicable precautions can be found:

- ▶ In the instructions for additional electrically operated devices (e.g. compressors for positioning systems, infusion pumps, nutrition probe ...)
- ▶ In the Standard DIN EN 60601-1-1:2002 (Safety of medical electrical systems)
- ▶ In the Standard VDE 0107:1994 (High power systems in hospitals)

■ The maximum long-term load for this bed insert is 250 kg (patient and accessories). Under this load, the bed insert can be operated without restriction.

■ The permitted weight of the patient depends on the total weight of the accessories attached at the same time (e.g. respirators, infusions,...)

#### Example:

Weight of accessories (with mattress):	Max. permitted weight of patient
10 kg	240 kg
40 kg	210 kg

■ This bed insert is not appropriate for patients less than 150 cm tall. Observe, especially with patients in poor health condition, the safety tips in chapter 4.

■ This bed insert may only be operated by trained personnel.

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■ This bed insert is useful for multiple reuse. Bear in mind the necessary requirements:

- ▶ See chapter 5;
- ▶ See chapter 6.

**This bed insert may only be operated under the operating instructions laid out in this instruction manual. Any other use is deemed to be inappropriate.**

**Sale of this product in North America, especially in the United States of America, is prohibited. Distribution and use of the bed insert, including distribution and use by third parties in these markets is forbidden by the Manufacturer.**

#### STRUCTURAL DESIGN

The bed insert is shipped disassembled. It can be transported without problems even in apartment homes. It consists of a reclining surface frame in two parts with mounted drive motors for

the head and foot sections. The reclining surface frame rests on a base frame. Two separate motors controls the up and down movements.

#### Reclining surface frame

The reclining surface frame has four sections: a movable backrest, a firm middle section, a moveable thigh rest and lower leg rest. The back- and thigh rests can be adjusted with the aid of electric motors. The horizontal height of the reclining surface can be adjusted. All adjustment functions are controlled by a hand switch.

#### Base frame

The base frame is equipped with height adjustable pedestal feet leveling on uneven floors. There are four fastening latches within the existing bed box. The base frame can be adapted to different bed widths.

#### Electric adjustment system

The electrical adjustment system of this bed insert is error protected, flame-retardant (V0) and consists of:

- an external transformer unit.
- The transformer unit consists of a power cable, the transformer, and a low voltage connection cable.
- The transformer creates a low voltage, which is safe for patient and user. The transformer supplies all drive motors with the safety low voltage through cables and a wiring harness. The connector on the base frame is moisture proof.
- the central control unit. All drive motors and the hand switch are hooked up to the central control unit by plug-in connections and work with the safety low-voltage;

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- a hand switch with a sturdy hook.
- The user can lock the adjustment functions of the hand switch if the poor clinical condition of the patient warrants it.
- two drive motors with integrated control for the horizontal height.
- one drive motor with integrated control for the thigh rest .
- one drive motor with integrated control for the backrest.

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#### DESCRIPTION OF MATERIALS

The bed insert is built for the most part from formed steel, the surface of which is covered with a polyester powder coating or has a metallic coating of zinc or chrome.

The spring wood slats consist of wood or wood materials, the surface which is sealed.

**All surfaces are harmless to skin contact.**

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### ASSEMBLY AND START UP PROCEDURE

The bed insert is shipped disassembled. Assembly must be done by two persons.

The bed insert ships in two cartons:

Carton 1

Divided reclining surface frame with motor drives, frame connectors, mattress brackets.

Carton 2

Base frame with motor, control box, hand switch and instruction manual.

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


#### ASSEMBLY

##### BASE FRAME AND RECLINING SURFACE FRAME


- ☞ Remove the contents of both cartons.
- ☞ Put the base frame [12] in the bed box. The connecting socket [14] must point to the head end of the bed box.
- ☞ Take the power adapter unit [15] out of its wrappings.
- ☞ Pull the plug out of the connecting socket [14].
- ☞ Plug the 4-poled plug (a) of the power adapter unit in this connecting socket. The 6-poled plug has a safety groove. Press the plug firmly into the connecting socket.
- ☞ Positioning the power adapter unit [15] secure in the bed box.

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☞ Plug the power plug of the power adapter  
into the wall socket.

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The two frame connectors [9] are secured to the  
top end reclining surface frame before leaving  
the factory.

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Please note, that the frame connectors for the  
right and left side are different.

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
☞ Remove the frame connectors from the  
reclining surface frame by loosening the knurled  
screws.

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☞ Insert the frame connectors into the tube  
ends at the top end.

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
The scissors of the base frame now drive up-  
wards as long as you press the button.

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☞ Move the scissors of the base frame to the  
top position.

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☞ Insert the halves of the reclining surface frame  
into each other.

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☞ Tighten the 4 knurled screws securely.  
Do not use pliers to tighten the screws!

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☞ Insert the reclining surface with the guide  
bars into the bearing rollers at the foot end.

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
☞ Push the reclining surface towards the foot end, until the frame connectors can be inserted into the location holes.

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☞ Secure the frame connectors each with a knurled screw.

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### CONNECTIONS OF MOTOR DRIVES

☞ Now disconnect the main power plug.

☞ Next connect the motor drives of the head and foot sections with the wiring harness [17].

☞ Insert the plug with the continuous cable into the connector for the motor drive [5] for the head section.

Secure the cable by snapping it tightly into the opening of the strain relief.

☞ Plug the last cable at the end of the wiring harness into the connector for the motor drive [3] of the foot section.

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Secure the cable by snapping it tightly into the opening of the strain relief.

☞ Now plug the power plug of the power adapter into the wall socket.

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#### Observe without fail:

When laying the cable in the base frame, make sure that through the up and down movement no shearing or squeezing damages the cables.

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

### MOUNTING OF MATTRESS BRACKETS

☞ Please screw the 5 mattress brackets to the reclining surface frames.

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## FITTING THE BED INSERT INTO EXISTING BEDS OR BED FRAMES

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☞ Check using the dimensional drawing whether the bed insert will fit into your bed or bed frame. To prevent any pinch points or shearing zones, please bear in mind that 25 mm safety clearance must be left all round between the bed and the bed insert.

☞ Lower the reclining surface frame all the way down in order to ascertain clearance within the bed box. Ensure there is free space equally on all sides.

☞ Raise the reclining surface frame into the highest position.

☞ Use the pedestal feet [13] to adjust to any unevenness of the floor. Turning the feet to the right will raise the frame, turning them to the left will lower it.

The bed insert can now be secured to the bed box.


☞ Loosen the 4 knurled screws near the fastening latches [11] and withdraw the fastening latches until they make contact with the bed box.

☞ Rotate the fastening latches up and use a wood screw to screw them to the bed box.

Choose the correct length of wood screws.

☞ Tighten the 2 knurled screws near the fastening latches securely! Do not use pliers to tighten the screws!

The bed insert is now firmly attached to the bed box.

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2,5 cm minimum tolerance  
between reclining surface  
frame and bed box

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The top of the reclining surface frame [1] must be 75 mm over the bed box!

The surrounding space between the bed insert and the bed box must be a minimum of 25 mm to avoid squeezing limbs while reaching into the bed box.

The mattress brackets [4] should not be below the top of the existing bed box. (Danger of injury!)

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## LOCATION REQUIREMENTS

- There must be enough room for the whole assemble of the bed insert. There must be no furniture, window sills, etc. in the way
- The drive area underneath the motor drives height adjustment must be kept free.
- A properly installed 230 volt wall plug must be available within the vicinity of the bed insert.
- Take care when attaching further accessory equipment (e.g compressors for positioning systems, etc) that their attachment is safe and their functionality is not impaired. Especially important is hereby also the safe placement of all moveable connecting cables, tubes, etc. If you have any further questions or need clarification contact the manufacturer of the accessory equipment or BURMEIER.

**Observe the following points when setting up the bed insert to minimize possible fire hazard through outside effects. Instruct the user and the authorized personnel in these points.**

- If possible use only fire-retardant mattresses and bed sheets.
- Only use appropriate firm mattresses according DIN 13014 with a density of 35 to 40 kg/m<sup>3</sup>. Further this mattress should be inflammable according to DIN 597 Section 1 and 2.
- Avoid smoking in bed because of the danger of fire.
- Use only technically faultless appliances (e.g. heating blankets) and other electronic devices ( e.g. lamps, radios).
- Take care that these devices are only used for the purpose intended and that they cannot accidentally get on top of or underneath the bed sheets (Fire hazard because of overheating).
- Absolutely avoid the use of extension cords or multi outlet extension cords under the bed (Leaking liquids can cause a fire hazard).
- Always plug the power plug directly into a wall socket. Extension cords and multi outlet extensions should not be used at all.

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**REACHING OPERATIONAL READINESS**

- After successful assembly complete an inspection according to chapter 6.
- Clean and disinfect the bed before first time use and before every further use according to chapter 5.

After successfully completing and heeding all steps on pages 12 to 17 the bed is now operationally ready.

**DISASSEMBLY**

- ☞ Raise the bed insert to the highest possible horizontal position.
- ☞ Remove the wood screws from the fastening latches of the bed box.
- ☞ Loosen the knurled screws on the fastening latches. Push the fastening latches all the way into the base frame. Secure the fastening latches with the knurled screws.
- ☞ Pull the power plug out of the wall socket.
- ☞ Pull out the plugs of the drive motors for the head and foot sections.
- ☞ Remove the knurled screws from the frame connectors.
- ☞ Lift up the reclining surface frame towards the head end.
- ☞ Pull the reclining surface frame towards the top end to remove them from the guide bars.
- ☞ Please put the reclining surface frame aside.
- ☞ Plug the power plug back into the wall socket.
- ☞ Lower the bed insert into its lowest horizontal position.
- ☞ Pull the plug out of the wall socket and remove the low voltage plug from the connector socket.
- ☞ Lift the base frame out of the bed box.
- ☞ If your transportation allows it leave the reclining frame assembled. Otherwise, separate the reclining surface frame.
- ☞ Put the frame connectors back into place on the reclining surface frame and screw them tightly.

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**SPECIAL SAFETY INSTRUCTIONS FOR THE DRIVE UNIT**

- Take care when adjusting the bed that no limbs of patients, users or any other persons, especially children, can be jammed or injured by getting under the movable parts or the bed frame.
- When using accessories with the bed it is necessary to: make sure that the arrangement of these accessories will not cause any places that might squeeze or hurt the patient. If this is not the case the operator must safely disable the adjustment function of the back- and thigh rests.
  - ▶ To do so place the hand switch out of the patient's reach (e.g. foot section) or lock the adjustment functions on the hand switch.
- Make sure that the power and hand switch cable cannot be jammed or damaged in other ways.
- Ensure that only hand switches that are in perfect condition are used. Faulty hand switches (e.g. with torn or indented foil) can result in spontaneous operation of the bed and can thus endanger the patient.

- Make sure that when adjusting the bed no obstacles such as furniture or slanted ceilings are in the way. This will prevent damage.
- Make sure that there are no objects on or in the spaces of the base frame.  
This could cause damage to the bed insert and/or the base frame.
- Take care when laying the hand switch cable that it is not squeezed by any pieces of the bed insert.
- Make sure that the connector socket is always closed with the plug when it is disconnected from the power adapter.
- Take care that the wiring harness or the cable of the hand switch is not pulled or squeezed in any way.

- ◆ An electronic overload detector switches off the central control unit, when the load is too great. After clearing the overload and pressing the hand switch again the central control unit will work again.

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- ◆ A continuous run time of 2 minutes should not be exceeded. Afterwards a pause of no less than 18 minutes should be observed. Another possibility: 1 minute continuous run time and 9 minutes pause.
- ◆ If the maximum continuous run time is not heeded (e.g. by constantly "playing" with the hand switch) a thermal safety component switches off the power for safety reasons. Depending on how large the over-load was, it can take several minutes before any more adjustments can be made.
- ◆ The adjustable range of every function is electrically/mechanically possible only within a predetermined range.
- ◆ Electrical disturbances affecting other nearby electronic appliances (e.g. crackling on the radio) cannot be excluded, despite adherence to all prescribed regulations. In such rare circumstances increase the distance between the appliances, do not use the same power socket or switch the affected appliance off for a while.

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## HAND SWITCH

The electrical adjustments can be made using the hand switch [8].

The adjustable range of every function is electrically/mechanically possible only within a predetermined range.

Locking mechanisms have been built into the hand switch for safety reasons. For the safety of the patient, adjustment functions can be locked at the hand switch, if, from a clinical point of view and with the approval of the attending physician, the condition of the patient warrants it.

- ◆ The hand switch, with its elastic hook can be hooked onto the bed.
- ◆ The spiral cable allows for greater mobility.
- ◆ The hand switch is water-tight and washable (IPX4).

### OPERATING THE HAND SWITCH

- Train the patient in the functions of the hand switch.
- The drive motors will run as long as the button is held down.
- The adjustment functions can be used in both directions.
- Basically the functions of the buttons are as shown :
  - raise
  - lower

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### Adjusting the backrest

Push this button if you want to adjust the backrest angle.

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### Adjusting the reclining surface height

Push this button if you want to adjust the height of the reclining surface.

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### Adjusting the thigh rest

Push this button to adjust the thigh rest angle.

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### Setting up the Auto-contour

Press this button if you want to adjust the angles of the back and thigh rest at the same time. This prevents the patient sliding to the foot end.

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### LOCKING FUNCTION

#### Only the operator is authorized to use the locking function!

In case the clinical condition of the patient is so critical that any adjustment of the beds by means of hand switch causes danger, the operator has to disable the bed immediately. The bed insert remains in the position it had at the time of disconnecting.

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Do not turn the locking key with force over the stop of the lock! The lock or the whole hand switch can be damaged.

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The locking key is delivered attached to the bed with a cable tie on the hand switch. The locking key is not intended for use by the patient. It must be removed from the hand switch. The operator or a person authorized by the physician should keep the locking key.

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
■ Turn the respective lock on the hand switch with the locking key clockwise in the locking position.

The color of the respective display will change from green to yellow.

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
Lock is vertical

Color of the display is: **green**

Buttons can be pushed ("Clicking sound")

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#### Drive disabled

Lock has been turned clockwise by about 15°

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## LOWER LEG REST

When the thigh rest is fixed, the lower leg rest [2] can be individually adjusted. For this reason there are two adjusting ratchets, which will lock the lower leg rest into position.

### Raising

☞ Raise the lower leg rest on the foot end - not with the mattress brackets - up to the desired level. The lower leg rest will automatically snap into place.

### Lowering

☞ First lift the lower leg rest all the way to the stopping point.  
☞ Slowly lower the lower leg rest.

There is a potential for injury if the lower leg rest falls straight down.

If you lower the thigh rest, the lower leg rest will be lowered automatically.

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## EMERGENCY LOWERING OF THE BACK REST

The upright back rest can be lowered manually in case of a power outage or failure of the power supply or the electrical drive systems.

**In this case two users are absolutely necessary!**

**Disregarding these safety guidelines and instructions may cause the backrest to fall out of control, which could lead to serious injuries for the user and the patient!**

Emergency lowering may only be carried out in the case of extreme emergency and by users who have a complete command of the procedures below.


We strongly advise you to practice emergency lowering of the backrest under normal conditions so that in case of an emergency you will be able to react quickly and properly.

☞ Before the emergency lowering, ease the load on the backrest.

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☞ The first user raised the back rest slightly by the outside edge of the head section and holds it in this position.

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- ☞ The second user now removes the spring pin from the lift rod. To do this swing away the bent bracket and pull the spring pin along with the bent bracket out of the lift rod of the back rest motor.
- ☞ The motor has now been detached from the back rest and swings down.
- ☞ Now the first user carefully lowers the back rest.

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**If the back rest falls out of control, the patient and the second user may be injured!**

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- ☞ To convert the bed back to its normal operating state: Swing up the connecting mount and secure the spring pin, then turn over the bent bracket.

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
### TRAPEZE POLE (ACCESSORIES)

- The maximum load of the trapeze pole is 75 kg at the front end.
- The trapeze pole is not suitable for rehabilitation training.

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The trapeze pole attached to the bed makes it easier for the patient to get into and out of bed.

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In both corners of the head section of the reclining surface frame are two trapeze pole sockets [8]. The trapeze pole socket has a notch at the top [c], which together with the metal pin [b] restricts the swivel range of the trapeze pole. The trapeze pole should be installed on the side where the patient gets in and out of bed.

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#### Attaching


- ☞ Insert the trapeze pole into the socket. The metal pin [b] must be in the notch [c] of the socket.

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#### Removal

- ☞ Pull the trapeze pole straight up out of the socket.

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#### Trapeze pole swivel range

Do not swivel the trapeze pole outside of the bed. There is a danger that through the pull, the bed will tip over.

The metal pin must always be in the notch of the socket.

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## EXTENDERS (ACCESSORIES)

In case the bed box requires the base frame to be higher, extenders for the pedestal feet are available.

Remove the feet from the pedestals, screw in the extenders and replace the feet under them.

## TRAPEZE

A trapeze (grab handle) can be attached to the trapeze pole, it is enclosed with the trapeze pole. The patient can hold on to the trapeze to make it easier to change to a different position.

Regularly check the trapeze at the strap for any damage (see Chapter 6). A damaged trapeze or strap must immediately be exchanged.

### Durability

There is a date clock on the trapeze. With normal use the trapeze will last for at least five years. Afterwards a sight check is to be performed periodically to determine if further use is permissible.

### Adjustable range of the trapeze

The height of the trapeze can be adjusted by means of the strap to accommodate a range of about 55 cm to 70 cm (measuring from the top of the mattress).

☞ Push the loop of the trapeze [14] over the first stud of the patient grab handle pole.



☞ Test by tightly pulling the trapeze downward that it is secure.

**Useful Tip:** The maximum lifting load of the trapeze pole is 75 kg at the front end.

■ The height of the trapeze can be adjusted by means of the strap.

■ Make sure the the strap is fed correctly through the buckle.

■ Make sure that the end of the strap is at least 3 cm above the buckle.

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


## BED EXTENSION

If an extra-long mattress is used, the bed insert can be equipped with a bed extension.

- ☞ Please remove the mattress brackets on the foot side.
- ☞ Please remove the caps on the foot side from the bed insert.

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- ☞ Secure the bed extension with the aid of screws.

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### GENERAL

Cleaning is the most important requirement for a successful chemical disinfection.

Generally a routine cleaning of the bed, if used by the same patient, is hygienically sufficient.

Disinfection of the bed frame is only necessary in case of visible contamination with infectious or potentially infectious material (blood, stool, pus) or by presence of an infectious disease under the direction of a physician.

**The bed insert must be cleaned and washed down with disinfectant in case of a patient change!**

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#### Before beginning the cleaning, please note:

- Unplug the power cable from the wall socket and store it in such a way that it does not come into excessive contact with water or other cleaning solutions (store in plastic bag).
- Make sure all plugs at the drive motors are properly connected.
- None of the electrical components must show any visible damage; otherwise water or cleaning solution could penetrate. This could lead to malfunction or damage to the electrical components.
- Before operating the bed again, make sure that there is no residual moisture on the electrical contacts, by drying or blowing on the power plug.
- The electrical components must not be cleaned with a water jet, a high pressure cleaner or other similar device!  
Clean only with a moist towel.
- If you suspect that water or any other form of moisture has penetrated the electrical components, pull the power plug out of the wall socket immediately. If already disconnected from the power supply, make sure it is not plugged in again. Clearly mark the bed insert as "Defective" and take it immediately out of service. Report this to the operator promptly.
- Failure to follow these directions could lead to major damage of the bed insert and its electrical components and lead to subsequent malfunctions.

## CLEANING AND DISINFECTION PLAN

- Remove the bed sheets and give them to the laundry.
- Clean all surfaces including the slatted boards and the reclining surface made of plastic inserts or wire mesh base with a mild and environmentally safe cleaning solution. This also applies to the hand switch.
- If the bed has visible contamination e.g. infectious or potentially infectious material has polluted it, it must be disinfected by wiping after cleaning. Use approved disinfectants for each surface according to the products listed by the DGHM (German Organization for Hygiene and Microbiology). This also applies to all beds from patients with illnesses subject to report according to §6 of the infection control law (German abbreviation: IfSG, Protection against Infection Act), colonization or infection with treatment resistant pathogens (e.g. MRSA, VRE) and all beds from the intensive care and infection

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control units. Use only the DGHM approved disinfectant in the approved concentration.

■ Disinfecting the pedestals is only necessary when visibly contaminated by infectious or potentially infectious material.

**Notice:** Continuous disinfection is only necessary for patients with a treatment resistant pathogen (e.g. MRSA) within the hospital.

## INSTRUCTING THE USER AND QUALIFIED PERSONNEL

To guarantee the correct procedure for cleaning and disinfection is followed, we recommend to correctly instruct the handler and qualified personnel. Thereby communicating that the following points need to be observed:

■ The clean bed insert is to be transported to the patient's home in such a way that it does not become polluted or contaminated.

■ When disassembling the bed insert we recommend cleaning and disinfecting it right away. The qualified personnel should be familiar with the steps necessary to thoroughly clean and disinfect (the operator will provide the individual steps necessary). It is necessary to use only the DGHM approved disinfectant in the approved concentration. The disinfectant must be appropriate for the surface being cleaned.

■ The qualified personnel must be equipped with water-tight one-time-use aprons and gloves.

■ Only fresh and clean towels must be used to clean and disinfect, these must then be given to the laundry.

■ After cleaning the bed insert, the qualified personnel must also disinfect their hands before going to the next patient. An appropriate disinfectant dispenser should be part of the qualified personnel's gear.

The immediate cleaning of the bed insert while disassembling has the benefit of not transporting "dirty" bed inserts or frame parts along with clean ones in the same truck. A transfer of potentially infectious spores which may be clinging to the used bed insert may be avoided. A transfer of spores of a hospitalacquired infection will be avoided with strict adherence to these recommendations.

If the bed inserts are to be stored for a while, they should be covered to keep them free of dust, accidental pollution and contamination.

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## CLEANING SUPPLIES AND DISINFECTANTS

Make use of these recommendations to keep the bed insert serviceable for a long time.

- The surfaces must be intact; any damage should be repaired immediately.
- We recommend a (damp) wash down. When choosing the cleaning solution make sure it is mild (gentle to skin and surfaces) and environmentally friendly. In general a regular household cleaner can be used.
- For cleaning the coated metal parts use a wet cloth with regular household cleaners.
- The disinfectant used for the wash down should be one approved for the purpose and listed by the DGHM (in the recommended dosage). Disinfectants and concentration listed by the Robert-Koch-Institute (the former Bundesgesundheitsamtes) should only be used in case of an epidemic according to the infection safety law and only on the orders of the appropriate health department.
- Despite the excellent mechanical resistance of the bed, scratches, markings, etc. which penetrate the coating should be resealed using a suitable medium to prevent the penetration of moisture. Contact BURMEIER or your specialized dealer with any questions.

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### Notice:

Surface disinfectants based on aldehyde generally have the advantage of working on a large array of pathogens, have a relatively low protein error and are environmentally friendly. Its main disadvantage is the allergy-causing and irritation potential. Compounds on the basis of glucoprotamine do not have this disadvantage and are very effective, though mostly more expensive. Disinfectants on the basis of chlorine splitting combinations, with prolonged contact or in too high doses, can be corrosive to metals, plastics, rubber and other materials. In addition to this, the compounds have a high protein error, are irritants to the mucous membranes and are not environmentally friendly.

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## CONTACT WITH DISINFECTANTS

- Follow the guidelines for each product and manufacturer. Observe the exact dosage! We recommend automatic dosage apparatus.
- Always mix the solution with cold water, to avoid vapors that aggravate the mucous membranes.
- Wear gloves to avoid direct skin contact.
- Do not store the mixed disinfectant in open containers with the cleaning rag inside. Close the container!
- Use closeable bottles with dosage pumps for wetting the cleaning rag.
- After finishing the surface disinfection ventilate the room.
- Use a wipe down disinfecting method, do not use a spray. Using the disinfectant as a spray causes the disinfectant to turn into a mist which could be inhaled.
- The wipe-down disinfection plays a major part.
- Do not use alcohol for larger areas.

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Operators of nursing beds are obligated to keep medical products in a safe state throughout their entire service life. This also includes regularly carrying out expert maintenance as well as safety checks.

If a bed insert is purchased for domestic use (not for commercial use), the regular safety inspection is not compulsory but is still recommended by the manufacturer.

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### Notice for Operators

This bed insert has been constructed and built to last safely for a long time. With appropriate operation and use, this bed has a life expectancy of 2 to 8 years. The life expectancy is based on the conditions and frequency of use.

Through repeated transport, assembly and disassembly, improper operation as well as long time use, it is possible that damage, defects and deterioration can occur. These deficiencies can lead to danger if they are not recognized and quickly fixed.

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Because of this there are legal bases for the regular enforcement of tests to ensure the safe condition of this medical product.

According to § 4 of the Medizinprodukte-Betreiberverordnung operators are obliged to maintain the product. Because of this the following regular inspection and function controls are to be performed by the operator as well as the user.

**Point out these necessary tests to the user.**

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### BY THE USER

Along with the regular substantial tests by technically qualified personnel, the normal user (nurses, supporting relatives, etc.) must also perform a sight and function check at shorter, regular intervals, as well as before a new patient.

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■ If damage or a malfunction is suspected, take the bed insert out of service immediately and disconnect it from the power until an exchange or repair can be made!

■ Contact your operator regarding exchange or repair of faulty parts.

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### Recommendation:

All electrical and mechanical components should be tested monthly. The main power cable as well as the hand switch cable should be checked after every mechanical load and after every relocation.

For this purpose use the checklist on the following **page 33**.

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**Checklist: Test protocol by the user**

Test	OK	Not OK	Defect description
<b>Visual check of electrical components</b>			
Hand switch and cable	Damage, Cable Routing		
Hand switch	Damage, Foil		
Power adapter unit	Damage, Cable Routing		
<b>Visual check of mechanical components</b>			
Trapeze pole/Socket	Damage, Deformation		
Bed frame	Damage, Deformation		
Reclining surface frame	Damage, Deformation		
<b>Function test of electrical components</b>			
Hand switch	Functiontest, Blocking function		
Strain relief, Main cable of power adapter	Tight fit		
<b>Function test of mechanical components</b>			
Emergency lowering of backrest	Test according to instruction manual		
Socket head cap screw	Tight fit		
Knurled screws	Tight fit		
Lower leg rest	Latching		
Accessories (e.g. trapeze pole, trapeze)	Secure, Damage		
<b>Signature of the tester:</b>	<b>Result of the test:</b>		<b>Date:</b>

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## BY THE OPERATOR

The operator of the bed insert according to MPBetreibV § 4 is obligated with every new assembly, every maintenance and throughout the operation, to perform regular inspections, to ensure the safe condition of the bed insert.

These tests are to be repeated during regular maintenance and depending on the use according to MPBetreibV § 4 and from the Directive of the German Employers Liability Insurance Association according to BGV A3.

- ◆ Observe this sequence for the inspection according to DIN EN 62353:
  - I. Visual inspection
  - II. Electric Measurements
  - III. Function Test
- ◆ The performance of the function test as well as the rating and documentation of the test results can, according to MPBetreibV § 4 only be performed by experts with the appropriate knowledge, training and experience as well as the tools necessary to perform it according to the rules.
- ◆ The performance of electrical measurements according to DIN EN 62353 can, with availability of proper measuring equipment, also be performed by an electro-technically trained person [according to BGV A3] with medical and appliance specific training.
- ◆ The assessment and documentation of the test results can only be performed by an expert with the corresponding knowledge, training and experience.

### Test method:

- Output power test: direct or differential according to DIN EN 62353
- Output power test according to the test unit manufacturers instructions.

### Limit:

- Output power less than 0.1mA.

### Test cycle:

We recommend a yearly visual inspection and performance check. If this checks was been ok, an electrical testing suffices every ten years.

If the 2% error quota (see BGV A3: §5, table 1 B) is shown to have been observed, the test cycle for inspections and functional testing can be extended to a maximum of two years.

If damage or a function error is suspected take the bed insert out of service immediately and disconnect it from the power until an exchange or repair can be made!

The test protocol draft on **pages 35 and 36** can be

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<b>Customer / med. Facility / practice:</b>						
<b>Address:</b>						
<b>Carried out:</b> <input type="checkbox"/> Repeat inspection			<input type="checkbox"/> Inspection prior initial operation			
<input type="checkbox"/> _____			<input type="checkbox"/> Inspection following repairs / servicing			
<b>Equipment type:</b> <input type="checkbox"/> Hospital bed <input type="checkbox"/> Care bed <input checked="" type="checkbox"/> Bed insert			<b>Protection classification:</b> <input type="checkbox"/> I <input checked="" type="checkbox"/> II			
<b>Bed type:</b> LIPPE 120, 140			<b>Inventory Number:</b>			
<b>Location:</b>			<b>Serial Number:</b>			
<b>Grade level (IEC60601-2-52):</b> <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5						
<b>Manufacturer:</b> Burmeier GmbH & Co. KG			<b>User-specific:</b> none			
<b>Testing equipment used (Type / Inventory Number):</b>			1.			
<b>MPG Classification:</b> Class I			2.			
<b>I. Visual Inspection</b>				<b>OK</b>	<b>Not OK</b>	<b>Defect Description</b>
<b>Visual inspection of the electrical components</b>						
Stickers and type plates	Presence, legible					
Housing of power adapter unit / Power adapter	Damage					
Motor housing and lifting rods of motors	Correct position, Damage					
Motor bolt	Check that the lock washer is properly in place					
Hand switch	Damage, Foil intact					
Motor-, hand switch- and main cables	Damage, Routing					
Plug and plug covering strip of control unit	Available, Correct position					
<b>Visual inspection of the mechanical components</b>						
Stickers and type plates	Presence, legible					
Safe working load / patient weight	Presence, legible					
Trapeze pole/socket	Damage, Deformation					
Bed frame	Damage, Deformation					
Anschlag mittels Aufnahmebolzen	Correct position					
Reclining surface	Damage, Deformation					
Welded seams	Split welded seams					
Knurled screws	Tight fit					
Wearing parts such as joints	Damage					
Locking ring of backrest/high rest	Check that the ring is properly in place					
<b>II. Electric measurements according to DIN EN 62353:2008-08 (leakage current, direct)</b>						
<b>Proceed as follow:</b>				<b>OK</b>	<b>Not OK</b>	<b>Defect Description</b>
1. Plug the power adapter into the test socket on the measuring instrument.						
2. Insert the plug of the 24 volt cable into the special adapter (extra accessory; request from BURMEIER).						
3. Connect the special adapter with the test probe of the test unit.						
4. Start procedures with the measuring instrument						
	Limit Value	Value of first measurement	Current actual value			
Result: Bed SK II (Type B)	0,1 mA	mA	mA			

III. Performance check		OK	Not OK	Defect Description
<b>Performance check of the electrical components</b>				
Motion limit switching of the motors	Automatic cut-off			
Hand switch, blocking functions	Test according to instructions for use			
Drive motors	Text, Abnormal Noise Development			
Strain relief and security of routing the cable	Secure, Correct position			
Strain relief for power adapter cable	Secure hold			
<b>Performance check of the mechanical components</b>				
Joints and pivots	Smooth operation			
Trapeze with belt strap	Secure hold			
Emergency lowering of the backrest	Test according to instructions for use			
Lower leg rest	Latching			
Accessories (e. g. trapeze pole, trapeze)	Secure, Damage, Suitability			
<b>Result of the inspection:</b>				
All values within permissible range: <input type="checkbox"/> yes <input type="checkbox"/> no		Next inspection date:		
Inspection passed: <input type="checkbox"/> yes <input type="checkbox"/> no				
<b>If inspection was not passed:</b> <input type="checkbox"/> Defective, do not use the bed! => Repair <input type="checkbox"/> Defective, do not use the bed! => Take out of Service <input type="checkbox"/> Bed does not meet the safety standards				
Test approval sticker applied: <input type="checkbox"/> yes <input type="checkbox"/> no				
Part of this inspection report:				
Comments:				
Inspected on:	Inspector:	Signature:		
Evaluated on:	Operator / Expert:	Signature:		

BURMEIER GmbH & Co KG		
Pivitsheider Straße 270 D - 32791 LAGE ☎ +49 - 05232 - 9841 - 0		
Model	LIPPE 120, 140	
Type	123456	
Serie	11 - 2002 - 123456	

Model	Bettenmodell
Type	Artikelnummer
Serie	Seriennummer

Model	Bed model
Type	Item Number
Serial	Serial Number

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## SPARE PARTS

The corresponding replacement parts can be obtained from BURMEIER by specifying the item number, the order number and the serial number. The necessary details are found on the type plate under the cross tubing of the reclining surface frame.

## MANUFACTURER'S ADDRESS

In order to maintain operational safety and the right to claim under warranty, only original BURMEIER replacement parts may be used!

For ordering replacement parts, customer service requests and other questions, please contact:

**Burmeier GmbH & Co. KG**  
(A subsidiary of the Stieglmeyer Group)

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## REPLACEMENT OF ELECTRICAL COMPONENTS

### Dangerous voltage, electric shock possible

- Before doing any work unplug the main power cable from the wall socket!
- Any work on the electrical equipment may only be performed by customer service, the motor manufacturer or qualified and authorized electricians under consideration of the applicable VDE rules and safety guidelines!

- Only dismantle the drive motors when the reclining surface is in the horizontal position, otherwise there is the risk of being injured by falling parts of the mattress base.

- All the drive components are maintenancefree and must not be opened. In the event of a malfunction, the corresponding components should always be replaced in full!

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### Exchange of the Hand switch

- ☞ Unplug the power plug from the wall socket.
- ☞ Remove the pull-out safety device on the control box.
- ☞ Unplug the plug at the end of the spiral cord of the hand switch from the control box connector.
- ☞ Lay the cable of the new hand switch so that it cannot be damaged by moving parts of the bed.
- ☞ Insert the plug at the end of the spiral cable of the new hand switch into the connector of the hand switch. Take care that the O-ring on the plug is not damaged, otherwise the connection may not be moisture-proof.
- ☞ Re-insert the pull-out safety device into the control box.
- ☞ Carry out a function test of the electrical adjustments after the exchange.

### Exchange of the power adapter

- ☞ To perform electrical tests the power adapter can be removed from the bed. The bed can remain occupied.
- ☞ Unplug the power cable from the wall socket.
- ☞ Unplug the plug of the 24-Volt connecting cable from the connecting socket at the head section.
- ☞ Connect the plug of the new power adapter into the connecting socket.
- ☞ Take care when laying the cable that it will not be damaged by moving parts or be run over.
- ☞ Carry out a function test of the electrical adjustments after the exchange.

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### TROUBLESHOOTING

The following table on page 40 is a guide to rectifying common malfunctions:

If a function fault is detected during operation, which the user cannot rectify with the use of the trouble-shooting table aid, the specialized staff for maintenance and repair of the appropriate owner is to be called.

The user may **under no circumstances** try to eliminate defects at the electrical equipment!

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### Dangerous voltage, electric shock possible!

Any work on the electrical equipment may only be performed by service personnel, the motor manufacturer or qualified and authorized electricians as required by the applicable VDE rules and safety guidelines!

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### Trouble-shooting table aid

Problem	Possible causes	Solution
Hand switch / drive system not functioning	<ul style="list-style-type: none"> <li>■ Hand switch plug not properly connected</li> <li>■ Transformer plug (24 volt) not properly connected</li> <li>■ Hand switch or drive system defective</li> <li>■ Functions are blocked on hand switch</li> </ul>	<ul style="list-style-type: none"> <li>☞ Test the plug connection</li> <li>☞ Plug the transformer plug into the connection socket</li> <li>☞ Inform your operator about necessary repairs</li> <li>☞ Unlock functions (see page 20)</li> </ul>
Green LED on power adapter not lit and drive systems are without function	<ul style="list-style-type: none"> <li>■ Power cable not properly connected</li> <li>■ Wall socket without power</li> <li>■ Power adapter defective</li> </ul>	<ul style="list-style-type: none"> <li>☞ Plug in the power cable properly</li> <li>☞ Test the wall socket and fuse box</li> <li>☞ Exchange the power adapter; inform your operator about necessary repairs</li> </ul>
Green LED on power adapter lit continuously, yellow LED not lit when any hand switch button is pushed	<ul style="list-style-type: none"> <li>■ Hand switch is defective</li> <li>■ Transformer plug (24 volt) not properly connected</li> <li>■ Functions are locked on hand switch</li> </ul>	<ul style="list-style-type: none"> <li>☞ Exchange the hand switch</li> <li>☞ Test the plug connection</li> <li>☞ Unlock the functions (see page 20)</li> </ul>
Drives operate in one direction only	<ul style="list-style-type: none"> <li>■ Hand switch or drive system is defective</li> </ul>	<ul style="list-style-type: none"> <li>☞ Inform your operator about necessary repairs</li> </ul>
Drives stop suddenly after prolonged adjustment times	<ul style="list-style-type: none"> <li>■ Heat switch in power adapter has been activated</li> </ul>	<ul style="list-style-type: none"> <li>☞ Too short minimum pause (see page 19)</li> <li>☞ Exchange the power adapter; inform your operator about necessary repairs</li> </ul>

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Only the Original BURMEIER accessories may be used. If other accessories are used, we can not be liable for any resulting injuries, defects and dangers.

When using side rail, drip/I.V., etc. with electrically adjustable beds, please note:

Make sure that the arrangement of accessories does not produce any crush or shearing zones for the patient when the back and leg rests are adjusted. If this cannot be ensured, the user must safely prevent the patient from adjusting the back- and leg rest. For this, lock the hand switch.

- ▶ To do so, attach the hand switch out of reach ( e.g. at the foot end of the bed), or lock hand switch adjustments.

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**MEASUREMENTS AND DIMENSIONS**
**Assembled bed insert:**

Total weight (size 200 x 90 cm)	
	Lippe 120 (200 x 120 cm)
	Lippe 140 (200 x 140 cm)
Safe load	
Patient weight, max.	

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**SPECIAL FEATURES**

- ◆ Electrical height adjustment of the reclining surface from about 29 to 74 cm.
- ◆ Electrical adjustment of the backrest from 0° to about 70°.
- ◆ Electrical adjustment of the leg rest from 0° to about 40°.
- ◆ Reclining surfaces in various sizes, divided into four areas (see page 5).
- ◆ Trapeze pole with trapeze handle available (accessories).

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### Motors for reclining surface height adjustment

Type
Power / Size / Lift
Input Voltage
Operating Time
Protection

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- ◆ Zulässiges Patientengewicht

### Explanation of Icons:

- ◆ Unit with Type B application according to IEC 601-1 (Special protection against electric shock)

- ◆ Protection Class II device, double insulation

- ◆ Transformer with thermal fuse block

- ◆ Safety transformer in accordance with VDE 05510

- ◆ Caution! Consult instruction manual

- ◆ The electrical equipment offers protection from internal dust deposits and is splash-proof from all sides

- ◆ The electrical equipment is dust-proof from internal dust deposits and is splash-proof from all sides

- ◆ Conformity mark in accordance with Medical Products directive 93/42 EWG

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### Noise level during motor adjustments

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The ambient conditions stated below must be observed:

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Storage temperature  
Relative humidity

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Ambient temperature  
Relative humidity  
Air pressure

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Risk analysis for medical devices
Safety for electromedical equipment
Electromagnetic compatibility
Medical electrical equipment: with fitness for use
Medical electrical equipment – Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Particular requirements for the safety and essential performance of medical beds
Medical electrical equipment: Particular requirements for the basic safety and essential performance of medical beds

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
hereby declare under sole responsibility as the manufacturer that the product model named below:

### Height-adjustable bed insert LIPPE 120, 140

in the version submitted complies with the regulations of the EC Directive 93/42/EEC for Medical Devices, last amended by Directive 2007/47/EC dated 5 September 2007. It is categorised as a Class I active medical device.

The relevant technical documentation is kept by the manufacturer's safety representative.

To evaluate the conformity to the Directives, all applicable parts of the following standards were referred to:

complies with the regulations of EC Guidelines 93/42/EEG for medical products and fulfills all additional safety requirements by the Highest State Authority from May 22, 2001 for care beds.


#### Harmonised standards:


- EN 14971 : 2013-04 Risk analysis for medical devices
- EN 60601-1 : 2007-07 Safety for electromedical equipment
- EN 60601-1-2 : 2007-12 Electromagnetic compatibility
- DIN EN 60601-1-6 : 2010-10 Medical electrical equipment: with fitness for use
- DIN EN 60601-1-11 : 2010-10 Medical electrical equipment – Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
- DIN EN 60601-2-52 : 2010-12 Particular requirements for the safety and essential performance of medical beds

#### International standards:

- IEC 60601-1-2-52 : 2009-12 Medical electrical equipment: Particular requirements for the basic safety and essential performance of medical beds

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No part of this manual may be reproduced with-  
out the prior written permission of the publisher!

All rights reserved.  
Subject to technical changes!

Stand: 2015.08.28

Updated: 2015.08.28